



eltB Antibody, Biotin conjugated

Product Code	CSB-PA315643LD01ENL
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0CK94
Immunogen	Recombinant Escherichia coli Heat-labile enterotoxin B chain protein (22-124AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Heat-labile enterotoxin B chain (LT-B, human) (LTH-B), eltB, ltpB
Immunogen Species	Escherichia coli
Research Area	Microbiology
Target Names	eltB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.