




mpt63 Antibody

Product Code	CSB-PA358712YA01MVZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WIP1
Immunogen	Recombinant Mycobacterium tuberculosis Immunogenic protein MPT63 protein (30-159AA)
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Immunogenic protein MPT63 (16 kDa immunoprotective extracellular protein) (Antigen MPT63), mpt63
Immunogen Species	Mycobacterium tuberculosis
Research Area	Others
Target Names	mpt63

Image



Western blot

All lanes: mpt63 antibody at 2μg/ml

Lane 1: Recombinant Mycobacterium tuberculosis Immunogenic protein MPT63 protein 1μg

Lane 2: Recombinant Mycobacterium tuberculosis Immunogenic protein MPT63 protein 0.1μg

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 17 kDa

Observed band size: 17 kDa