



# mpt63 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA358712YD01MVZ
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	P9WIP1
<b>Immunogen</b>	Recombinant Mycobacterium tuberculosis Immunogenic protein MPT63 protein (30-159AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Mycobacterium tuberculosis
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Immunogenic protein MPT63 (16 kDa immunoprotective extracellular protein) (Antigen MPT63), mpt63
<b>Immunogen Species</b>	Mycobacterium tuberculosis (strain ATCC 25618 / H37Rv)
<b>Research Area</b>	Others
<b>Target Names</b>	mpt63