



mpt63 Antibody, Biotin conjugated

Product Code	CSB-PA358712YD01MVZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P9WIP1
Immunogen	Recombinant Mycobacterium tuberculosis Immunogenic protein MPT63 protein (30-159AA)
Raised In	Rabbit
Species Reactivity	Mycobacterium tuberculosis
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Immunogenic protein MPT63 (16 kDa immunoprotective extracellular protein) (Antigen MPT63), mpt63
Immunogen Species	Mycobacterium tuberculosis (strain ATCC 25618 / H37Rv)
Research Area	Others
Target Names	mpt63
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.