




PDZD11 Antibody

Product Code	CSB-PA707840HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5EBL8
Immunogen	Recombinant Human PDZ domain-containing protein 11 protein (1-140AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	PDZ domain-containing protein 11 (ATPase-interacting PDZ protein) (Plasma membrane calcium ATPase-interacting single-PDZ protein) (PMCA-interacting single-PDZ protein), PDZD11, AIPP1 PDZK11 PISP
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PDZD11

Image



Immunohistochemistry of paraffin-embedded human appendix tissue using CSB-PA707840HA01HU at dilution of 1:100