



CSF2 Antibody, HRP conjugated

Product Code	CSB-PA06475B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04141
Immunogen	Recombinant Human Granulocyte-macrophage colony-stimulating factor protein (18-144AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Granulocyte-macrophage colony-stimulating factor (GM-CSF) (Colony-stimulating factor) (CSF) (Molgramostin) (Sargramostim), CSF2, GMCSF
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CSF2