



Timp1 Antibody, HRP conjugated

Product Code	CSB-PA16719B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30120
Immunogen	Recombinant Rat Metalloproteinase inhibitor 1 protein (24-217AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Metalloproteinase inhibitor 1 (Tissue inhibitor of metalloproteinases 1) (TIMP-1), Timp1, Timp-1
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Cardiovascular
Target Names	Timp1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.