



PIP Antibody, HRP conjugated

Product Code	CSB-PA018020LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12273
Immunogen	Recombinant Human Prolactin-inducible protein (29-146AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Prolactin-inducible protein (Gross cystic disease fluid protein 15) (GCDFP-15) (Prolactin-induced protein) (Secretory actin-binding protein) (SABP) (gp17), PIP, GCDFP15 GPI4
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	PIP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.