



gD Antibody, HRP conjugated

Product Code	CSB-PA18929B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P57083
Immunogen	Recombinant Human herpesvirus 1 Envelope glycoprotein D protein (26-339AA)
Raised In	Rabbit
Species Reactivity	Human herpesvirus 1
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Envelope glycoprotein D (gD), gD, US6
Immunogen Species	Human herpesvirus 1
Research Area	Others
Target Names	gD
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.