



CXCL8 Antibody, Biotin conjugated

Product Code	CSB-PA011671YD01DO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P41324
Immunogen	Recombinant Dog Interleukin-8 protein (23-101AA)
Raised In	Rabbit
Species Reactivity	Dog
Tested Applications	ELISA
Relevance	IL-8 is a chemotactic factor that attracts neutrophils, basophils, and T-cells, but not monocytes. It is also involved in neutrophil activation. It is released from several cell types in response to an inflammatory stimulus.
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Interleukin-8 (IL-8) (C-X-C motif chemokine 8) (Chemokine (C-X-C motif) ligand 8), CXCL8, IL8
Species	Dog
Research Area	Cardiovascular
Target Names	CXCL8