




GAGE5 Antibody

Product Code	CSB-PA009186EA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13069
Immunogen	Recombinant Human G antigen 5 protein (1-117AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	G antigen 5 (GAGE-5) (Cancer/testis antigen 4.5) (CT4.5), GAGE5
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	GAGE5

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA009186EA01HU at dilution of 1:100