



FGF1 Antibody

Product Code	CSB-PA06379A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05230
Immunogen	Recombinant Human Fibroblast growth factor 1 protein (16-155AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fibroblast growth factor 1 (FGF-1) (Acidic fibroblast growth factor) (aFGF) (Endothelial cell growth factor) (ECGF) (Heparin-binding growth factor 1) (HBGF-1), FGF1, FGFA
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	FGF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.