



CCL16 Antibody

Product Code	CSB-PA08967A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O15467
Immunogen	Recombinant Human C-C motif chemokine 16 protein (24-120AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	C-C motif chemokine 16 (Chemokine CC-4) (HCC-4) (Chemokine LEC) (IL-10-inducible chemokine) (LCC-1) (Liver-expressed chemokine) (Lymphocyte and monocyte chemoattractant) (LMC) (Monotactin-1) (MTN-1) (NCC-4) (Small-inducible cytokine A16), CCL16, ILINCK NCC4 SCYA16
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CCL16
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.