



crp Antibody, HRP conjugated

Product Code	CSB-PA359785LB01ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0ACJ8
Immunogen	Recombinant Escherichia coli cAMP-activated global transcriptional regulator CRP protein (1-210AA)
Raised In	Rabbit
Species Reactivity	Escherichia coli
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	cAMP-activated global transcriptional regulator CRP (Catabolite activator protein) (CAP) (Catabolite gene activator) (cAMP receptor protein) (CRP) (cAMP regulatory protein), crp, cap csm
Immunogen Species	Escherichia coli (strain K12)
Research Area	Others
Target Names	crp
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.