



S100a9 Antibody, Biotin conjugated

Product Code	CSB-PA020642HD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P31725
Immunogen	Recombinant Mouse Protein S100-A9 protein (2-113AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein S100-A9 (Calgranulin-B) (Leukocyte L1 complex heavy chain) (Migration inhibitory factor-related protein 14) (MRP-14) (p14) (S100 calcium-binding protein A9), S100a9, Cagb Mrp14
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	S100a9
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.