



TRMT112 Antibody, HRP conjugated

Product Code	CSB-PA883416LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UI30
Immunogen	Recombinant Human tRNA methyltransferase 112 homolog protein (1-125AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Multifunctional methyltransferase subunit TRM112-like protein (tRNA methyltransferase 112 homolog), TRMT112
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TRMT112
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.