



TRMT112 Antibody, FITC conjugated

Product Code	CSB-PA883416LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UI30
Immunogen	Recombinant Human tRNA methyltransferase 112 homolog protein (1-125AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Multifunctional methyltransferase subunit TRM112-like protein (tRNA methyltransferase 112 homolog), TRMT112
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TRMT112