



BMP4 Antibody, HRP conjugated

Product Code	CSB-PA05799B0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P12644
Immunogen	Recombinant Human Bone morphogenetic protein 4 protein (303-408AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Bone morphogenetic protein 4 (BMP-4) (Bone morphogenetic protein 2B) (BMP-2B), BMP4, BMP2B DVR4
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	BMP4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.