



ipaH9.8 Antibody, HRP conjugated

| | |
|----------------------------|--|
| Product Code | CSB-PA513237HB01SZC |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | D2AJU0 |
| Immunogen | Recombinant Shigella flexneri serotype X E3 ubiquitin-protein ligase ipaH9.8 protein (1-545AA) |
| Raised In | Rabbit |
| Species Reactivity | Shigella flexneri serotype X |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | E3 ubiquitin-protein ligase ipaH9.8 (EC 2.3.2.27) (Invasion plasmid antigen ipaH9.8) (RING-type E3 ubiquitin transferase ipaH9.8), ipaH9.8 |
| Immunogen Species | Shigella flexneri serotype X (strain 2002017) |
| Research Area | Others |
| Target Names | ipaH9.8 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |