



ipaH9.8 Antibody, Biotin conjugated

Product Code	CSB-PA513237HD01SZC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	D2AJU0
Immunogen	Recombinant Shigella flexneri serotype X E3 ubiquitin-protein ligase ipaH9.8 protein (1-545AA)
Raised In	Rabbit
Species Reactivity	Shigella flexneri serotype X
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase ipaH9.8 (EC 2.3.2.27) (Invasion plasmid antigen ipaH9.8) (RING-type E3 ubiquitin transferase ipaH9.8), ipaH9.8
Immunogen Species	Shigella flexneri serotype X (strain 2002017)
Research Area	Others
Target Names	ipaH9.8
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.