




Cd82 Antibody

Product Code	CSB-PA004961LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P40237
Immunogen	Recombinant Mouse CD82 antigen protein (111-227AA)
Raised In	Rabbit
Species Reactivity	Mouse, Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CD82 antigen (C33 antigen) (IA4) (Inducible membrane protein R2) (Metastasis suppressor Kangai-1 homolog) (CD antigen CD82), Cd82, Kai1
Immunogen Species	Mus musculus (Mouse)
Research Area	Immunology
Target Names	Cd82

Image



Western blot
All lanes: Cd82 antibody at 2µg/ml
Lane 1: Jurkat whole cell lysate
Lane 2: A549 whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 30 kDa
Observed band size: 30 kDa