




NR1H4 Antibody

Product Code	CSB-PA839406LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96R11
Immunogen	Recombinant Human Bile acid receptor protein (1-476AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Purification Method	Antigen affinity purification
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	NR1H4

Image



Western Blot

Positive WB detected in: HEK293 whole cell lysate

All lanes: NR1H4 antibody at 3µg/ml
Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 56, 55, 49 kDa

Observed band size: 56 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.