




CALM1 Antibody

Product Code	CSB-PA08285A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P62158
Immunogen	Recombinant Human Calmodulin protein (2-149AA)
Raised In	Rabbit
Species Reactivity	Human, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CaM, CALM1, CALM, CAM, CAM1, CALM2, CAM2, CAMB, CALM3, CALML2, CAM3, CAMC, CAMIII
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	CALM1

Image



Western blot
 All lanes: CALM1 antibody at 2µg/ml + Mouse skeletal muscle tissue
 Secondary
 Goat polyclonal to rabbit IgG at 1/10000 dilution
 Predicted band size: 17 kDa
 Observed band size: 17 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.