



HA Antibody, Biotin conjugated

Product Code	CSB-PA18939D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	R4NN21
Immunogen	Recombinant Influenza A virus Hemagglutinin protein (19-524AA)
Raised In	Rabbit
Species Reactivity	Influenza A virus(H7N9)
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Hemagglutinin [Cleaved into: Hemagglutinin HA1 chain; Hemagglutinin HA2 chain], HA
Immunogen Species	Influenza A virus(H7N9)
Research Area	Microbiology
Target Names	HA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.