



ompA Antibody, Biotin conjugated

Product Code	CSB-PA16479D0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8Z7S0
Immunogen	Recombinant Salmonella typhi Outer membrane protein A protein (27-349AA)
Raised In	Rabbit
Species Reactivity	Salmonella typhi
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Outer membrane protein A, ompA
Immunogen Species	Salmonella typhi
Research Area	Others
Target Names	ompA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.