



Cma1 Antibody

Product Code	CSB-PA005599LA01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21844
Immunogen	Recombinant Mouse Chymase protein (22-246AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Relevance	Major secreted protease of mast cells with suspected roles in vasoactive peptide generation, extracellular matrix degradation, and regulation of gland secretion.
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Chymase (EC 3.4.21.39) (Alpha-chymase) (Mast cell chymase 1) (Mast cell protease 5) (mMCP-5) (Mast cell protease I), Cma1, Mcpt5
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Cma1