



Cma1 Antibody, Biotin conjugated

Product Code	CSB-PA005599LD01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21844
Immunogen	Recombinant Mouse Chymase protein (22-246AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Chymase (EC 3.4.21.39) (Alpha-chymase) (Mast cell chymase 1) (Mast cell protease 5) (mMCP-5) (Mast cell protease I), Cma1, Mcpt5
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Cma1