




VP3 Antibody

Product Code	CSB-PA860325DA01CID
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q99152
Immunogen	Recombinant Chicken anemia virus Apoptin protein (1-121AA)
Raised In	Rabbit
Species Reactivity	Chicken anemia virus
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Apoptin, VP3
Immunogen Species	Chicken anemia virus
Research Area	Others
Target Names	VP3

Image



Western blot

All lanes: Chicken anemia virus apoptin antibody at 2μg/ml + GST-Apoptin fusion protein at 0.4μg
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 40 kDa
Observed band size: 40 kDa