



cry1Ab Antibody

Product Code	CSB-PA363358HA01BDC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0A370
Immunogen	Recombinant Bacillus thuringiensis subsp. kurstaki Pesticidal crystal protein Cry1Ab protein (1022-1155AA)
Raised In	Rabbit
Species Reactivity	Bacillus thuringiensis subsp. kurstaki
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pesticidal crystal protein Cry1Ab (130 kDa crystal protein) (Crystalline entomocidal protoxin) (Insecticidal delta-endotoxin CryIA(b)), cry1Ab, bt2 cry-1-2 cry1A(b) cryIA(b) cryIC1
Immunogen Species	Bacillus thuringiensis subsp. kurstaki
Research Area	Others
Target Names	cry1Ab
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.