



MBNL1 Antibody, HRP conjugated

Product Code	CSB-PA889088LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NR56
Immunogen	Recombinant Human Muscleblind-like protein 1 protein (1-382AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Muscleblind-like protein 1 (Triplet-expansion RNA-binding protein), MBNL1, EXP KIAA0428 MBNL
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	MBNL1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.