



LTB Antibody, HRP conjugated

Product Code	CSB-PA771253LB01EQV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q862Z7
Immunogen	Recombinant Chimpanzee Lymphotoxin-beta protein (49-244AA)
Raised In	Rabbit
Species Reactivity	Chimpanzee
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lymphotoxin-beta (LT-beta) (Tumor necrosis factor C) (TNF-C) (Tumor necrosis factor ligand superfamily member 3), LTB, TNFC TNFSF3
Immunogen Species	Pan troglodytes (Chimpanzee)
Research Area	Others
Target Names	LTB