



LTB Antibody, HRP conjugated

Product Code	CSB-PA771253LB01EQV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q862Z7
Immunogen	Recombinant Chimpanzee Lymphotoxin-beta protein (49-244AA)
Raised In	Rabbit
Species Reactivity	Chimpanzee
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lymphotoxin-beta (LT-beta) (Tumor necrosis factor C) (TNF-C) (Tumor necrosis factor ligand superfamily member 3), LTB, TNFC TNFSF3
Immunogen Species	Pan troglodytes (Chimpanzee)
Research Area	Others
Target Names	LTB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.