



Msln Antibody, HRP conjugated

Product Code	CSB-PA015044LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q61468
Immunogen	Recombinant Mouse Mesothelin protein (298-600AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mesothelin (Pre-pro-megakaryocyte-potentiating factor) [Cleaved into: Megakaryocyte-potentiating factor (MPF); Mesothelin, cleaved form], Msln, Mes Mpf
Immunogen Species	Mus musculus (Mouse)
Research Area	Tags & Cell Markers
Target Names	Msln
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.