




Cd74 Antibody

Product Code	CSB-PA004956LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P10247
Immunogen	Recombinant Rat H-2 class II histocompatibility antigen gamma chain protein (57-280AA)
Raised In	Rabbit
Species Reactivity	Rat, Mouse
Tested Applications	ELISA, WB; Recommended dilution: WB:1:1000-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	H-2 class II histocompatibility antigen gamma chain (Ia antigen-associated invariant chain) (Ii) (MHC class II-associated invariant chain) (CD antigen CD74), Cd74
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Immunology
Target Names	Cd74

Image



Western blot

All lanes: Cd74 antibody at 2µg/ml + Mouse lung tissue

Secondary

Goat polyclonal to rabbit IgG at 1/10000 dilution

Predicted band size: 32, 25 kDa

Observed band size: 32 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.