





COX4I1 Antibody

Product Code	CSB-PA00509A0Rb
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P13073
Immunogen	Recombinant Human Cytochrome c oxidase subunit 4 isoform 1, mitochondrial protein (23-169AA)
Raised In	Rabbit
Species Reactivity	Human, Rat, Zebrafish
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Cytochrome c oxidase subunit 4 isoform 1, mitochondrial (Cytochrome c oxidase polypeptide IV) (Cytochrome c oxidase subunit IV isoform 1) (COX IV-1), COX4I1, COX4
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	COX4I1

Image



Western blot
All lanes: COX4I1 antibody at 10µg/ml
Lane 1: Rat muscle tissue
Lane 2: Zebrafish lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 20 kDa
Observed band size: 20 kDa



Immunohistochemistry of paraffin-embedded human gastric cancer using CSB-PA00509A0Rb at dilution of 1:100