




TMX3 Antibody

Product Code	CSB-PA846645LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96JJ7
Immunogen	Recombinant Human Protein disulfide-isomerase TMX3 protein (25-195AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein disulfide-isomerase TMX3 (EC 5.3.4.1) (Thioredoxin domain-containing protein 10) (Thioredoxin-related transmembrane protein 3), TMX3, KIAA1830 TXNDC10
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TMX3

Image



Immunohistochemistry of paraffin-embedded human testis tissue using CSB-PA846645LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.