



MRPL58 Antibody, HRP conjugated

Product Code	CSB-PA613490HB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14197
Immunogen	Recombinant Human Peptidyl-tRNA hydrolase ICT1, mitochondrial protein (30-206AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Peptidyl-tRNA hydrolase ICT1, mitochondrial (EC 3.1.1.29) (39S ribosomal protein L58, mitochondrial) (MRP-L58) (Digestion subtraction 1) (DS-1) (Immature colon carcinoma transcript 1 protein) (Mitochondrial large ribosomal subunit protein mL62), MRPL58, DS1 ICT1
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	MRPL58
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.