




GMPR Antibody

Product Code	CSB-PA339756HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P36959
Immunogen	Recombinant Human GMP reductase 1 protein (1-345AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GMP reductase 1 (GMPR 1) (EC 1.7.1.7) (Guanosine 5'-monophosphate oxidoreductase 1) (Guanosine monophosphate reductase 1), GMPR, GMPR1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	GMPR

Image



Immunohistochemistry of paraffin-embedded human liver cancer using CSB-PA339756HA01HU at dilution of 1:100

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.