




PLD5 Antibody

Product Code	CSB-PA836686LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N7P1
Immunogen	Recombinant Human Inactive phospholipase D5 protein (92-536AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, IHC; Recommended dilution: IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Inactive phospholipase D5 (Inactive PLD 5) (Inactive choline phosphatase 5) (Inactive phosphatidylcholine-hydrolyzing phospholipase D5) (PLDc), PLD5
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PLD5

Image



Immunohistochemistry of paraffin-embedded human breast cancer using CSB-PA836686LA01HU at dilution of 1:100

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.