



ERGIC3 Antibody, HRP conjugated

Product Code	CSB-PA896688LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y282
Immunogen	Recombinant Human Endoplasmic reticulum-Golgi intermediate compartment protein 3 protein (47-341AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Endoplasmic reticulum-Golgi intermediate compartment protein 3 (Serologically defined breast cancer antigen NY-BR-84), ERGIC3, C20orf47 ERV46 SDBCAG84
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ERGIC3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.