



PTPN6 Antibody, FITC conjugated

Product Code	CSB-PA019043HC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29350
Immunogen	Recombinant Human Tyrosine-protein phosphatase non-receptor type 6 protein (1-597AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Tyrosine-protein phosphatase non-receptor type 6 (EC 3.1.3.48) (Hematopoietic cell protein-tyrosine phosphatase) (Protein-tyrosine phosphatase 1C) (PTP-1C) (Protein-tyrosine phosphatase SHP-1) (SH-PTP1), PTPN6, HCP PTP1C
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PTPN6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.