





STAMPB Antibody

Product Code	CSB-PA022794HA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O95630
Immunogen	Recombinant Human STAM-binding protein (1-424AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB, IHC; Recommended dilution: WB:1:1000-1:5000, IHC:1:20-1:200
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	STAM-binding protein (EC 3.4.19.-) (Associated molecule with the SH3 domain of STAM) (Endosome-associated ubiquitin isopeptidase), STAMPB, AMSH
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	STAMPB


Image



Western blot
All lanes: STAMPB antibody at 12µg/ml
Lane 1: Jurkat whole cell lysate
Lane 2: Hela whole cell lysate
Secondary
Goat polyclonal to rabbit IgG at 1/10000 dilution
Predicted band size: 49, 40 kDa
Observed band size: 49 kDa



Immunohistochemistry of paraffin-embedded human prostate cancer using CSB-PA022794HA01HU at dilution of 1:100



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using CSB-PA022794HA01HU at dilution of 1:100