



E6 Antibody, HRP conjugated

Product Code	CSB-PA365854LB01HML
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03126
Immunogen	Recombinant Human papillomavirus type 16 Protein E6 protein (1-158AA)
Raised In	Rabbit
Species Reactivity	Human papillomavirus type 16
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein E6, E6
Immunogen Species	Human papillomavirus type 16
Research Area	Microbiology
Target Names	E6