



E6 Antibody, Biotin conjugated

Product Code	CSB-PA365854LD01HML
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03126
Immunogen	Recombinant Human papillomavirus type 16 Protein E6 protein (1-158AA)
Raised In	Rabbit
Species Reactivity	Human papillomavirus type 16
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein E6, E6
Immunogen Species	Human papillomavirus type 16
Research Area	Microbiology
Target Names	E6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.