



Slc39a13 Antibody

Product Code	CSB-PA644804LA01RA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2M1K6
Immunogen	Recombinant Rat Zinc transporter ZIP13 protein (130-233AA)
Raised In	Rabbit
Species Reactivity	Rat
Tested Applications	ELISA
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc transporter ZIP13 (Solute carrier family 39 member 13) (Zrt- and Irt-like protein 13) (ZIP-13), Slc39a13, Zip13
Immunogen Species	Rattus norvegicus (Rat)
Research Area	Others
Target Names	Slc39a13