



Slc39a13 Antibody

| | |
|----------------------------|----------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-PA644804LA01RA |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q2M1K6 |
| Immunogen | Recombinant Rat Zinc transporter ZIP13 protein (130-233AA) |
| Raised In | Rabbit |
| Species Reactivity | Rat |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Zinc transporter ZIP13 (Solute carrier family 39 member 13) (Zrt- and Irt-like protein 13) (ZIP-13), Slc39a13, Zip13 |
| Immunogen Species | Rattus norvegicus (Rat) |
| Research Area | Others |
| Target Names | Slc39a13 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |