



KLRG1 Antibody, HRP conjugated

Product Code	CSB-PA846610LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96E93
Immunogen	Recombinant Human Killer cell lectin-like receptor subfamily G member 1 protein (60-195AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Killer cell lectin-like receptor subfamily G member 1 (C-type lectin domain family 15 member A) (ITIM-containing receptor MAFA-L) (MAFA-like receptor) (Mast cell function-associated antigen), KLRG1, CLEC15A MAFA MAFAL
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	KLRG1