



CD7 Antibody, FITC conjugated

Product Code	CSB-PA004953YC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P09564
Immunogen	Recombinant Human T-cell antigen CD7 protein (26-180AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	T-cell antigen CD7 (GP40) (T-cell leukemia antigen) (T-cell surface antigen Leu-9) (TP41) (CD antigen CD7), CD7
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.