



CA12 Antibody, HRP conjugated

Product Code	CSB-PA004367YB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O43570
Immunogen	Recombinant Human Carbonic anhydrase 12 protein (25-301AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Carbonic anhydrase 12 (EC 4.2.1.1) (Carbonate dehydratase XII) (Carbonic anhydrase XII) (CA-XII) (Tumor antigen HOM-RCC-3.1.3), CA12
Immunogen Species	Homo sapiens (Human)
Research Area	Cancer
Target Names	CA12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.