



BID Antibody, FITC conjugated

Product Code	CSB-PA002698HC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P55957
Immunogen	Recombinant Human BH3-interacting domain death agonist protein (1-195AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	BH3-interacting domain death agonist (p22 BID) (BID) [Cleaved into: BH3-interacting domain death agonist p15 (p15 BID); BH3-interacting domain death agonist p13 (p13 BID); BH3-interacting domain death agonist p11 (p11 BID)], BID
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	BID
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.